

المرحلة الرابعة / طفيليات عملي (٢)

- **Protozoa:** unicellular organisms, e.g. Plasmodium (malaria)
- **Metazoa:** multicellular organisms, e.g. helminths (worms) and arthropods (ticks, lice)

Diagnosis of Parasitic Infections

1. Clinical
2. Laboratory

Purpose of laboratory diagnosis :

- Confirmation of clinical suspicion
- Identification of unsuspected infection

Specimens:


- ❖ Stool
- ❖ Blood
- ❖ Serum and plasma
- ❖ Others (anal swab, duodenal aspirate, sputum, urine, urogenital specimen)
- ❖ Tissues.

Stool examination:

Sample collection:

- Sample is collected in clean, dry container
- Handled carefully
- Sometimes use preservative (10% formalin)
- Samples in some cases fresh (amoeba)
- Liquid and soft stool examined within 15 min
- Not mixed with urine or disinfectant (as they will kill trophozoites)

المرحلة الرابعة/ طفيليات عملي (٢)



Morphology of *Dientamoeba fragilis* from a stool sample. Trophozoites exhibit an ameba-like morphology and are often bi-nucleated.

- Specimens obtained by enema or laxatives are often positive for worm eggs or adult worm.

Examination of the stool sample:

Gross examination:

- Mucooid blood stained (acute amoebic dysentery), Parasites can be detected (nematodes, cestodes)

Microscopic examination:

- Saline mount
- Iodine Mount
- Thick smears – not commonly used
- **Permanent stained smears**
 - Iron hematoxylene
 - Whearley's trichrome stain
- Concentration methods
 - Floatation techniques
 - Sedimentation techniques

Antigen detection

Molecular diagnosis

Microscopic examination

Direct wet mount:


- Thin emulsion of small amount of faeces
- **Few drops of saline**

المرحلة الرابعة / طفيليات عملي (٢)


- Sometimes add lugol's iodine (nuclear details, glycogen vacuole in cyst)
- Protozoa (trophozoite), cyst, eggs and larva of helminths, crystals

Concentration methods:

- Scanty parasites in the sample
- **Floatation** (eggs and cyst float , solution of high specific gravity)
 1. saturated sodium chloride (ascaris, hookworms)
 2. Zinc sulphate centrifugation floatation (cyst, nematodes).
- **Sedimentation** (solution of low specific gravity): formol ether



Formol Ether Sedimentation



- Ether adsorbs fecal debris & floats.
- Formalin fixes & preserves the specimen.